

MAJnMAR
MEDICAL PRODUCTS

Products Catalogue

MAJnMAR Medical Products

MAJnMAR Medical Products · 203-1130 Austin Ave, Coquitlam, BC, V3K 3P5, CANADA **www.majnmar.com**

MAJnMAR MEDICAL PRODUCTS
203-1130 Austin Ave
Coquitlam, British Columbia
V3K 3P5, CANADA
Tel: +1 778 872 3876
Fax: +1 604 939 9944

CONTENTS

Index

01 Trocars 3

Bladeless 4
Optical 5
Autoshield 4

02 Endoscopic Linear Cutting Staplers & Reloads 6

Optimum Stapler 6
Optimal Stapler 8
Optimum Reloads 7
Optimal Reloads 9

03 Powered Staplers & Reloads 10

Powered Endoscopic Linear Cutting Stapler & Reloads 10
E-Force I Reloads 11
Reloads/Optimum 13
Charger/ Battery cover for E-Force II powered stapler 15
Optimal Reloads 16
E-Force I Stapler 11
E-Force II Staplers 12
Optimal Reloads 14
E-Force III 15

04 Hernia Staplers & Suturing 17

PremiaTack Endoscopic Hernia Stapler 17
EndoKnit Suturing Device 18
PremiaTack A Absorbable Endoscopic Hernia Stapler 17

05 Clip Appliers & Clips 19

OptiClip 19
PremiaClip 20
PolyClip 21
Reusable Titanium Clip Appliers 22
PolyClip Applier 19
Reusable Polymer Clip Applier 20
Reusable Clip Remover 21
Titanium Clips 22

06 Laparoscopic Instruments & Accessories 23

Veress Needle 23
Endobag 23
Endoscopic Instruments 24
Coagulation Electrode 26
Suction irrigation 23
Needle Holder 24
Disposable Instruments 25
Blunt Tip Dissector, Knot Pusher, Camera Sleeve 26

07 Open Surgery Staplers 27

Linear Cutting Staplers & Reloads 27
Circular Stapler 29
Linear Staplers & Reloads 28
Hemorrhoids Stapler 30

08 Silicone Drainage System 31

Round channeled Drain 32
Flat Channeled Drain 33
Silicone Bulb reservoir 34
Silicone T-Tube Drain 35
Round Perforated Drain 33
Flat Perforated Drain 34
I -Vac Reservoir 35
Calibration Tube 36

MAJnMAR MEDICAL PRODUCTS
203-1130 Austin Ave
Coquitlam, British Columbia
V3K 3P5, CANADA
Tel: +1 778 872 3876
Fax: +1 604 939 9944

PRODUCT RANGE

Products at a Glance



Trocars

Page 3



Optimum Stapler & Reloads

Page 6



Optimal Stapler & Reloads

Page 8



E-Force I Stapler & Reloads

Page 11



E-Force II Staplers

Page 12



E-Force III

Page 15



PremiaTack Hernia Staplers

Page 17



EndoKnit Suturing Device

Page 18



OptiClip Clip Appliers

Page 19



PolyClip & PremiaClip Appliers

Page 20



PolyClip Polymer Clips

Page 21



Reusable Clip Remover & Appliers

Page 21



Veress Needle

Page 23



Suction Irrigation

Page 23



Endobag

Page 23




Needle Holder

Page 24



Endoscopic Instruments

Page 24



Coagulation Electrode

Page 26



Linear Cutting Staplers & Reloads

Page 27



Linear Staplers & Reloads

Page 28



Circular Stapler

Page 29



Hemorrhoids Stapler

Page 30



Silicone Drains

Page 31



Silicone Bulb & I-Vac Reservoirs

Page 34




Calibration Tube

Page 36

MAJnMAR MEDICAL PRODUCTS
203-1130 Austin Ave
Coquitlam, British Columbia
V3K 3P5, CANADA
Tel: +1 778 872 3876
Fax: +1 604 939 9944

Trocars



MAJnMAR produces high quality Trocars of all types, Bladeless, Auto Shield and Optical available in multiple sizes (3, 5, 10, 12 and 15mm).

FEATURES:

- Insufflation optimally maintained by double sealers.
- Stability threads provide easier insertion with maximum fixation.
- Clear Cannula for increased vision and safety with beveled tip for easier extraction of small specimen.
- Ergonomic Cannula design to provide solid grip and easy insertion with curved smooth edges to prevent Patients or Surgeons injuries.

MAJnMAR MEDICAL PRODUCTS
203-1130 Austin Ave
Coquitlam, British Columbia
V3K 3P5, CANADA
Tel: +1 778 872 3876
Fax: +1 604 939 9944

Bladeless



FEATURES:

- Blunt Obturator tip with dual tissue separators to minimize the fascial defect through muscle-splitting and to allow surgeons to forego closure of fascia.

Type	Ref#	Description	Packaging
Set	MNM-303B	i-Port Disposable Bladeless Trocar Set, 3mm/65mm	6pcs/Box
	MNM-505BS	i-Port Disposable Bladeless Trocar Set, 5mm/70mm	6pcs/Box
	MNM-505B	i-Port Disposable Bladeless Trocar Set, 5mm/110mm	6pcs/Box
	MNM-510B	i-Port Disposable Bladeless Trocar Set, 5-10mm/110mm	6pcs/Box
	MNM-512B	i-Port Disposable Bladeless Trocar Set, 5-12mm/110mm	6pcs/Box
	MNM-505BL	i-Port Disposable Bladeless Trocar Set Long, 5mm/150mm	6pcs/Box
	MNM-510BL	i-Port Disposable Bladeless Trocar Set Long, 5-10mm/150mm	6pcs/Box
	MNM-512BL	i-Port Disposable Bladeless Trocar Set, Long, 5-12mm/150mm	6pcs/Box
Kit	MNM-505BSK	i-Port Disposable Bladeless Trocar Kit, 1 Obturator 2 cannulas, 5mm/70mm	10pcs/Box
	MNM-505Bk	i-Port Disposable Bladeless Trocar Kit, 1 Obturator 2 Cannulas, 5mm/110mm	10pcs/Box
	MNM-510Bk	i-Port Disposable Bladeless Trocar Kit, 1 Obturator 2 Cannulas, 5-10mm/110mm	10pcs/Box
	MNM-512Bk	i-Port Disposable Bladeless Trocar Kit, 1 Obturator 2 Cannulas, 5-12mm/110mm	10pcs/Box

Autoshield



FEATURES:

- Flat Bladed Obturator provides a smooth and safe insertion and features.
- Auto retracting blade shield which prevents damage of internal structures during trocar insertion.

Type	Ref#	Description	Packaging
Set	MNM-505S	i-Port Disposable Autoshield Trocar Set, 5mm/110mm	6pcs/Box
	MNM-510S	i-Port Disposable Autoshield Trocar Set, 5-10mm/110mm	6pcs/Box
	MNM-512S	i-Port Disposable Autoshield Trocar Set, 5-12mm/110mm	6pcs/Box
Kit	MNM-505SK	i-Port Disposable Autoshield Trocar Kit, 1 Obturator 2 Cannulas, 5mm/110mm	10pcs/Box
	MNM-510SK	i-Port Disposable Autoshield Trocar Kit, 1 Obturator 2 Cannulas, 5-10mm/110mm	10pcs/Box
	MNM-512SK	i-Port Disposable Autoshield Trocar Kit, 1 Obturator 2 Cannulas, 5-12mm/110mm	10pcs/Box

MAJnMAR MEDICAL PRODUCTS
203-1130 Austin Ave
Coquitlam, British Columbia
V3K 3P5, CANADA
Tel: +1 778 872 3876
Fax: +1 604 939 9944

Optical



FEATURES:

- Clear Dolphin tip Obturator offers outstanding visualization during Trocar insertion; provided with easily operated Scope locking system.

Type	Ref#	Description	Packaging
Set	MNM-505V	i-Port Disposable Optical Trocar Set, 5mm/110mm	6pcs/Box
	MNM-510V	i-Port Disposable Optical Trocar Set, 5-10mm/110mm	6pcs/Box
	MNM-512V	i-Port Disposable Optical Trocar Set, 5-12mm/110mm	6pcs/Box
	MNM-515V	i-Port Disposable Optical Trocar Set, 5-15mm/110mm	6pcs/Box
Kit	MNM-505VK	i-Port Disposable Optical Trocar Kit, 1 Obturator 2 Cannulas, 5mm/110mm	10pcs/Box
	MNM-510VK	i-Port Disposable Optical Trocar Kit, 1 Obturator 2 Cannulas, 5-10mm/110mm	10pcs/Box
	MNM-512VK	i-Port Disposable Optical Trocar Kit, 1 Obturator 2 Cannulas, 5-12mm/110mm	10pcs/Box
Set	MNM-505VH	i-Port Disposable Optical Trocar with handle Set, 5-12mm/110mm	6pcs/Box
	MNM-510VH	i-Port Disposable Optical Trocar with handle Set, 5-12mm/110mm	6pcs/Box
	MNM-512VH	i-Port Disposable Optical Trocar with handle Set, 5-12mm/110mm	6pcs/Box

MAJnMAR MEDICAL PRODUCTS
203-1130 Austin Ave
Coquitlam, British Columbia
V3K 3P5, CANADA
Tel: +1 778 872 3876
Fax: +1 604 939 9944

Endoscopic Linear Cutting Staplers & Reloads



MAJnMAR Endoscopic Linear Cutting Staplers & Reloads enable surgeons to operate with confidence while handling the broadest range of tissues and surgical applications. We are committed to delivering outstanding stapling performance that assures patient safety and surgeon's peace of mind.

FEATURES:

- Compression Force: 30 N.
- Burst Pressure: 5.0 kPa vs. normal Intra-Abdominal Pressure (IAP) of 0.67–0.93 kPa, with a maximum of 2 kPa. This results in a safety factor of 2.5 to 7 times.
- New Knife with Each Reload: ensures a precise, clean cut, reducing the risk of disease transmission in cancer procedures.
- Knife and Staples Mechanism: advances 15mm with each firing, requiring only 4 firings to complete the stapling process of a 60mm reload, resulting in saving time and improving surgical efficiency.
- Smooth Firing Action: provides a smoother firing experience compared to competing devices.

Optimum Stapler

FEATURES:

- Ergonomic Design: Compact and lightweight for optimal maneuverability and ease of use in confined spaces.
- Safety Mechanism: Advanced features to enhance patient safety and reduce the risk of complications.
- Multiple Reload Options: Supports a range of reload sizes, offering flexibility for various surgical needs.



Ref#	Size	Description	Packaging
MNM-ECL	260	Long Shaft, 260mm length, 12mm	1 PC/BOX
MNM-ECM	160	Medium Shaft, 160mm length, 12mm	1 PC/BOX
MNM-ECS	60	Short Shaft, 60mm length, 12mm	1 PC/BOX

MAJnMAR MEDICAL PRODUCTS
203-1130 Austin Ave
Coquitlam, British Columbia
V3K 3P5, CANADA
Tel: +1 778 872 3876
Fax: +1 604 939 9944

Optimum Reloads

FEATURES:

- Provides 3 staggered staple lines for a secure and reliable anastomosis.
- Delivers unified staples height for all staples.
- Features color-coded packaging for easy product identification and quick access.



Ref#	Size	Number of staples	Staples Height	Closure Height	Color	Packaging
MNM-ECR3020A	30	48	2	0.7	Articulating, 30mm, Grey	6 PCS/BOX
MNM-ECR3025A	30	48	2.5	1.0	Articulating, 30mm, White	6 PCS/BOX
MNM-ECR3035A	30	48	3.5	1.5	Articulating, 30mm, Blue	6 PCS/BOX
MNM-ECR4520A	45	66	2	0.7	Articulating, 45mm, Grey	6 PCS/BOX
MNM-ECR4525A	45	66	2.5	1.0	Articulating, 45mm,White	6 PCS/BOX
MNM-ECR4535A	45	66	3.5	1.5	Articulating, 45mm, Blue	6 PCS/BOX
MNM-ECR4540A	45	66	4	1.8	Articulating, 45mm, Gold	6 PCS/BOX
MNM-ECR4548A	45	66	4.8	2.0	Articulating, 45mm, Green	6 PCS/BOX
MNM-ECR6020A	60	90	2	0.7	Articulating, 60mm, Grey	6 PCS/BOX
MNM-ECR6025A	60	90	2.5	1.0	Articulating, 60mm, White	6 PCS/BOX
MNM-ECR6035A	60	90	3.5	1.5	Articulating, 60mm, Blue	6 PCS/BOX
MNM-ECR6040A	60	90	4	1.8	Articulating, 60mm, Gold	6 PCS/BOX
MNM-ECR6048A	60	90	4.8	2.0	Articulating, 60mm, Green	6 PCS/BOX
MNM-ECR3020	30	48	2	0.7	Straight, 30mm, Grey	6 PCS/BOX
MNM-ECR3025	30	48	2.5	1.0	Straight, 30mm, White	6 PCS/BOX
MNM-ECR3035	30	48	3.5	1.5	Straight, 30mm, Blue	6 PCS/BOX
MNM-ECR4520	45	66	2	0.7	Straight, 45mm, Grey	6 PCS/BOX
MNM-ECR4525	45	66	2.5	1.0	Straight, 45mm,White	6 PCS/BOX
MNM-ECR4535	45	66	3.5	1.5	Straight, 45mm, Blue	6 PCS/BOX
MNM-ECR4540	45	66	4	1.8	Straight, 45mm, Gold	6 PCS/BOX
MNM-ECR4548	45	66	4.8	2.0	Straight, 45mm, Green	6 PCS/BOX
MNM-ECR6020	60	90	2	0.7	Straight, 60mm, Grey	6 PCS/BOX
MNM-ECR6025	60	90	2.5	1.0	Straight, 60mm, White	6 PCS/BOX
MNM-ECR6035	60	90	3.5	1.5	Straight, 60mm, Blue	6 PCS/BOX
MNM-ECR6040	60	90	4	1.8	Straight, 60mm, Gold	6 PCS/BOX
MNM-ECR6048	60	90	4.8	2.0	Straight, 60mm, Green	6 PCS/BOX

MAJnMAR MEDICAL PRODUCTS
203-1130 Austin Ave
Coquitlam, British Columbia
V3K 3P5, CANADA
Tel: +1 778 872 3876
Fax: +1 604 939 9944

Optimal Stapler

FEATURES:



- Clean Cut Line: Provides a precise and smooth cut, minimizing tissue trauma.
- Ergonomic Design: Compact and lightweight for optimal maneuverability and ease of use in confined spaces.
- Safety Mechanism: Advanced features to enhance patient safety and reduce the risk of complications.
- Multiple Reload Options: Supports a range of reload sizes, offering flexibility for various surgical needs.


Ref#	Size	Descripion	Packaging
MNM-ECLT	260	Long Shaft, 260mm length, 12mm	1 PC/BOX
MNM-ECMT	160	Medium Shaft, 160mm length, 12mm	1 PC/BOX
MNM-ECST	60	Short Shaft, 60mm length, 12mm	1 PC/BOX

MAJnMAR MEDICAL PRODUCTS
203-1130 Austin Ave
Coquitlam, British Columbia
V3K 3P5, CANADA
Tel: +1 778 872 3876
Fax: +1 604 939 9944

Optimal Reloads

FEATURES:

- Color-Coded System: Utilizes a color-coded system for easy identification of staple height and reload type, streamlining the selection process during surgery.
- Provides 3 staggered staple lines for a secure and reliable anastomosis.
- Incorporates a gradient of staple heights, from larger to smaller, moving from outer to inner rows, optimizing blood flow to the cut line and facilitating faster healing.




Ref#	Size	Number of Staples	Staples Height	Description	Packaging
MNM-ECR30GA	30	48	a=1.0 b= 1.5 c= 2.0	Articulating, 30mm, Grey	6 PCS/BOX
MNM-ECR30TA	30	48	a=2.0 b= 2.5 c= 3.0	Articulating, 30mm, Tan	6 PCS/BOX
MNM-ECR30PA	30	48	a=3.0 b= 3.5 c= 4.0	Articulating, 30mm, Purple	6 PCS/BOX
MNM-ECR45GA	45	66	a=1.0 b= 1.5 c= 2.0	Articulating, 45mm, Grey	6 PCS/BOX
MNM-ECR45TA	45	66	a=2.0 b= 2.5 c= 3.0	Articulating, 45mm, Tan	6 PCS/BOX
MNM-ECR45PA	45	66	a=3.0 b= 3.5 c= 4.0	Articulating, 45mm, Purple	6 PCS/BOX
MNM-ECR45BA	45	66	a=4.0 b= 4.5 c= 5.0	Articulating, 45mm, Black	6 PCS/BOX
MNM-ECR60GA	60	90	a=1.0 b= 1.5 c= 2.0	Articulating,60mm, Grey	6 PCS/BOX
MNM-ECR60TA	60	90	a=2.0 b= 2.5 c= 3.0	Articulating,60mm, Tan	6 PCS/BOX
MNM-ECR60PA	60	90	a=3.0 b= 3.5 c= 4.0	Articulating,60mm, Purple	6 PCS/BOX
MNM-ECR60BA	60	90	a=4.0 b= 4.5 c= 5.0	Articulating,60mm, Black	6 PCS/BOX
MNM-ECR30GAT	30	48	a=1.0 b= 1.5 c= 2.0	Articulating, 30mm, Grey, Tilted	6 PCS/BOX
MNM-ECR30TAT	30	48	a=2.0 b= 2.5 c= 3.0	Articulating, 30mm, Tan, Tilted	6 PCS/BOX
MNM-ECR30PAT	30	48	a=3.0 b= 3.5 c= 4.0	Articulating, 30mm, Purple, Tilted	6 PCS/BOX
MNM-ECR45GAT	45	66	a=1.0 b= 1.5 c= 2.0	Articulating, 45mm, Grey, Tilted	6 PCS/BOX
MNM-ECR45TAT	45	66	a=2.0 b= 2.5 c= 3.0	Articulating, 45mm, Tan, Tilted	6 PCS/BOX
MNM-ECR45PAT	45	66	a=3.0 b= 3.5 c= 4.0	Articulating, 45mm, Purple, Tilted	6 PCS/BOX
MNM-ECR45BAT	45	66	a=4.0 b= 4.5 c= 5.0	Articulating, 45mm, Black, Tilted	6 PCS/BOX
MNM-ECR60GAT	60	90	a=1.0 b= 1.5 c= 2.0	Articulating,60mm, Grey, Tilted	6 PCS/BOX
MNM-ECR60TAT	60	90	a=2.0 b= 2.5 c= 3.0	Articulating,60mm, Tan, Tilted	6 PCS/BOX
MNM-ECR60PAT	60	90	a=3.0 b= 3.5 c= 4.0	Articulating,60mm, Purple, Tilted	6 PCS/BOX
MNM-ECR60BAT	60	90	a=4.0 b= 4.5 c= 5.0	Articulating,60mm, Black, Tilted	6 PCS/BOX

MAJnMAR MEDICAL PRODUCTS
203-1130 Austin Ave
Coquitlam, British Columbia
V3K 3P5, CANADA
Tel: +1 778 872 3876
Fax: +1 604 939 9944

Powered Endoscopic Linear Cutting Stapler & Reloads







Powered Staplers provide Surgeons with peace of mind and ensures patient safety through its solid yet smooth firing mechanism, resulting in perfect staple formation and clean cutline, hence ideal anastomosis.

MAJnMAR MEDICAL PRODUCTS
203-1130 Austin Ave
Coquitlam, British Columbia
V3K 3P5, CANADA
Tel: +1 778 872 3876
Fax: +1 604 939 9944

E-Force I Stapler



Ref#	Description	Packaging
MNM-EF13445X	Length 340MM, Reload size 45 MM, Mechanical Articulation	1 PC/BOX
MNM-EF13445Y	Length 340MM, Reload size 45 MM, Adaptive Articulation	1 PC/BOX
MNM-EF13460X	Length 340MM, Reload size 60 MM, Mechanical Articulation	1 PC/BOX
MNM-EF13460Y	Length 340MM, Reload size 60 MM, Adaptive Articulation	1 PC/BOX
MNM-EF14445X	Length 440MM, Reload size 45 MM, Mechanical Articulation	1 PC/BOX
MNM-EF14445Y	Length 440MM, Reload size 45 MM, Adaptive Articulation	1 PC/BOX
MNM-EF14460X	Length 440MM, Reload size 60 MM, Mechanical Articulation	1 PC/BOX
MNM-EF14460Y	Length 440MM, Reload size 60 MM, Adaptive Articulation	1 PC/BOX

E-Force I Reloads

FEATURES:



- Color-coded packaging for easy product identification and quick.
- Ease of Loading: Designed for quick and easy loading into the E-FORCE I Powered Stapler, minimizing downtime during surgery.
- Offers an affordable surgical stapling solution that strikes the perfect

balance between performance and affordability.

Ref#	Number of staples	Closure Heights	Description	Packaging
MNM-ECR45RP	64	0.7	Length 45mm, Staple height 2.0mm, Gray	12 PCS/BOX
MNM-ECR45WP	64	1.0	Length 45mm, Staple height 2.6mm,White	12 PCS/BOX
MNM-ECR45BP	64	1.5	Length 45mm,Staple height 3.5mm, Blue	12 PCS/BOX
MNM-ECR45YP	64	1.8	Length 45mm, Staple height 3.8mm, Gold	12 PCS/BOX
MNM-ECR45GP	64	2.0	Length 45mm, Staple height 4.1mm, Green	12 PCS/BOX
MNM-ECR45KP	64	2.3	Length 45mm, Staple height 4.3mm, Black	12 PCS/BOX
MNM-ECR60RP	88	0.7	Length 60mm, Staple height 2.0mm, Gray	12 PCS/BOX
MNM-ECR60WP	88	1.0	Length 60mm, Staple height 2.6mm, White	12 PCS/BOX
MNM-ECR60BP	88	1.5	Length 60mm, Staple height 3.5mm, Blue	12 PCS/BOX
MNM-ECR60YP	88	1.8	Length 60mm, Staple height 3.8mm, Gold	12 PCS/BOX
MNM-ECR60GP	88	2.0	Length 60mm, Staple height 4.1mm, Green	12 PCS/BOX
MNM-ECR60KP	88	2.3	Length 60mm, Staple height 4.3mm, Black	12 PCS/BOX

MAJnMAR MEDICAL PRODUCTS
203-1130 Austin Ave
Coquitlam, British Columbia
V3K 3P5, CANADA
Tel: +1 778 872 3876
Fax: +1 604 939 9944

E-Force II Staplers



FEATURES:

- Equipped with a digital chip scan screen to verify the authenticity of each reload.
- Customizable Options: Available in various sizes and configurations to meet the specific needs of different surgical procedures.
- Efficient Stapling Speed: Delivers rapid and consistent staple deployment.
- Available with a rechargeable battery option to support environmental sustainability.

Ref#	Description	Packaging
MNM-EF2L	E-FORCE 2 Large	1 PC/BOX
MNM-EF2M	E-FORCE 2 Medium	1 PC/BOX
MNM-EF2S	E-FORCE 2 Small	1 PC/BOX
MNM-EF2L-R	E-FORCE 2 Large with Rechargeable Battery	1 PC/BOX
MNM-EF2M-R	E-FORCE 2 Medium with Rechargeable Battery	1 PC/BOX
MNM-EF2S-R	E-FORCE 2 Small with Rechargeable Battery	1 PC/BOX

MAJnMAR MEDICAL PRODUCTS
203-1130 Austin Ave
Coquitlam, British Columbia
V3K 3P5, CANADA
Tel: +1 778 872 3876
Fax: +1 604 939 9944


Reloads/Optimum



Ref#	Size	Number of Staples	Staples Height	Closure Height	Description	Packaging
MNM-ECR3020AC	30	48	2	0.7	Articulating, 30mm, Grey with Chip	6 PCS/BOX
MNM-ECR3025AC	30	48	2.5	1.0	Articulating, 30mm, White with Chip	6 PCS/BOX
MNM-ECR3035AC	30	48	3.5	1.5	Articulating, 30mm, Blue with Chip	6 PCS/BOX
MNM-ECR4520AC	45	66	2	0.7	Articulating, 45mm, Grey with Chip	6 PCS/BOX
MNM-ECR4525AC	45	66	2.5	1.0	Articulating, 45mm, White with Chip	6 PCS/BOX
MNM-ECR4535AC	45	66	3.5	1.5	Articulating, 45mm, Blue with Chip	6 PCS/BOX
MNM-ECR4540AC	45	66	4	1.8	Articulating, 45mm, Gold with Chip	6 PCS/BOX
MNM-ECR4548AC	45	66	4.8	2.0	Articulating, 45mm, Green with Chip	6 PCS/BOX
MNM-ECR6020AC	60	90	2	0.7	Articulating, 60mm, Grey with Chip	6 PCS/BOX
MNM-ECR6025AC	60	90	2.5	1.0	Articulating, 60mm, White with Chip	6 PCS/BOX
MNM-ECR6035AC	60	90	3.5	1.5	Articulating, 60mm, Blue with Chip	6 PCS/BOX
MNM-ECR6040AC	60	90	4	1.8	Articulating, 60mm, Gold with Chip	6 PCS/BOX
MNM-ECR6048AC	60	90	4.8	2.0	Articulating, 60mm, Green with Chip	6 PCS/BOX
MNM-ECR3020C	30	48	2	0.7	Straight, 60mm, Grey with Chip	6 PCS/BOX
MNM-ECR3025C	30	48	2.5	1.0	Straight, 60mm, White with Chip	6 PCS/BOX
MNM-ECR3035C	30	48	3.5	1.5	Straight, 60mm, Blue with Chip	6 PCS/BOX
MNM-ECR4520C	45	66	2	0.7	Straight, 60mm, Grey with Chip	6 PCS/BOX
MNM-ECR4525C	45	66	2.5	1.0	Straight, 60mm, White with Chip	6 PCS/BOX
MNM-ECR4535C	45	66	3.5	1.5	Straight, 60mm, Blue with Chip	6 PCS/BOX
MNM-ECR4540C	45	66	4	1.8	Straight, 60mm, Gold with Chip	6 PCS/BOX
MNM-ECR4548C	45	66	4.8	2.0	Straight, 60mm, Green with Chip	6 PCS/BOX
MNM-ECR6020C	60	90	2	0.7	Straight, 60mm, Grey with Chip	6 PCS/BOX
MNM-ECR6025C	60	90	2.5	1.0	Straight, 60mm, White with Chip	6 PCS/BOX
MNM-ECR6035C	60	90	3.5	1.5	Straight, 60mm, Blue with Chip	6 PCS/BOX
MNM-ECR6040C	60	90	4	1.8	Straight, 60mm, Gold with Chip	6 PCS/BOX
MNM-ECR6048C	60	90	4.8	2.0	Straight, 60mm, Green with Chip	6 PCS/BOX

MAJnMAR MEDICAL PRODUCTS
203-1130 Austin Ave
Coquitlam, British Columbia
V3K 3P5, CANADA
Tel: +1 778 872 3876
Fax: +1 604 939 9944

Optimal Reloads



Ref #	Size	Number of staples	Staples Height	Description	Packaging
MNM-ECR30GAC	30	48	a=1.0 b= 1.5 c= 2.0	Articulating, 30mm, Grey, with Chip	6 PCS/BOX
MNM-ECR30TAC	30	48	a=2.0 b= 2.5 c= 3.0	Articulating, 30mm, Tan, with Chip	6 PCS/BOX
MNM-ECR30PAC	30	48	a=3.0 b= 3.5 c= 4.0	Articulating, 30mm, Purple, with Chip	6 PCS/BOX
MNM-ECR45GAC	45	66	a=1.0 b= 1.5 c= 2.0	Articulating, 45mm, Grey, with Chip	6 PCS/BOX
MNM-ECR45TAC	45	66	a=2.0 b= 2.5 c= 3.0	Articulating, 45mm, Tan, with Chip	6 PCS/BOX
MNM-ECR45PAC	45	66	a=3.0 b= 3.5 c= 4.0	Articulating, 45mm, Purple, with Chip	6 PCS/BOX
MNM-ECR45BAC	45	66	a=4.0 b= 4.5 c= 5.0	Articulating, 45mm, Black, with Chip	6 PCS/BOX
MNM-ECR60GAC	60	90	a=1.0 b= 1.5 c= 2.0	Articulating, 60mm, Grey, with Chip	6 PCS/BOX
MNM-ECR60TAC	60	90	a=2.0 b= 2.5 c= 3.0	Articulating, 60mm, Tan, with Chip	6 PCS/BOX
MNM-ECR60PAC	60	90	a=3.0 b= 3.5 c= 4.0	Articulating, 60mm, Purple, with Chip	6 PCS/BOX
MNM-ECR60BAC	60	90	a=4.0 b= 4.5 c= 5.0	Articulating, 60mm, Black, with Chip	6 PCS/BOX
MNM-ECR30GATC	30	48	a=1.0 b= 1.5 c= 2.0	Articulating, 30mm, Grey, Tilted Tip, with Chip	6 PCS/BOX
MNM-ECR30TATC	30	48	a=2.0 b= 2.5 c= 3.0	Articulating, 30mm, Tan, Tilted Tip, with Chip	6 PCS/BOX
MNM-ECR30PATC	30	48	a=3.0 b= 3.5 c= 4.0	Articulating, 30mm, Purple, Tilted Tip, with Chip	6 PCS/BOX
MNM-ECR45GATC	45	66	a=1.0 b= 1.5 c= 2.0	Articulating, 45mm, Grey, Tilted Tip, with Chip	6 PCS/BOX
MNM-ECR45TATC	45	66	a=2.0 b= 2.5 c= 3.0	Articulating, 45mm, Tan, Tilted Tip, with Chip	6 PCS/BOX
MNM-ECR45PATC	45	66	a=3.0 b= 3.5 c= 4.0	Articulating, 45mm, Purple, Tilted Tip, with Chip	6 PCS/BOX
MNM-ECR45BATC	45	66	a=4.0 b= 4.5 c= 5.0	Articulating, 45mm, Black, Tilted Tip, with Chip	6 PCS/BOX
MNM-ECR60GATC	60	90	a=1.0 b= 1.5 c= 2.0	Articulating, 60mm, Grey, Tilted Tip, with Chip	6 PCS/BOX
MNM-ECR60TATC	60	90	a=2.0 b= 2.5 c= 3.0	Articulating, 60mm, Tan, Tilted Tip, with Chip	6 PCS/BOX
MNM-ECR60PATC	60	90	a=3.0 b= 3.5 c= 4.0	Articulating, 60mm, Purple, Tilted Tip, with Chip	6 PCS/BOX
MNM-ECR60BATC	60	90	a=4.0 b= 4.5 c= 5.0	Articulating, 60mm, Black, Tilted Tip, with Chip	6 PCS/BOX

MAJnMAR MEDICAL PRODUCTS
203-1130 Austin Ave
Coquitlam, British Columbia
V3K 3P5, CANADA
Tel: +1 778 872 3876
Fax: +1 604 939 9944

Charger/ Battery cover for E-Force II powered stapler

Battery charger



- Main voltage: 5V
- Protection
- Charging indicator, in the battery
- Battery is not included
- Charging currency: 1A

Battery cover

- Engineered to shield a rechargable battery to maintain a sterile field.



Ref#	Descripion	Packaging
MNM-BB1	charger for E-FORCE II	1 PC/BOX
MNM-BC	battery cover E-Force II	1 PC/BOX

E-Force III

FEATURES:




- Equipped with a digital chip scan screen to verify the authenticity of each reload.
- Customizable Options: Available in various sizes and configurations to meet the specific needs of different surgical procedures.
- Efficient Stapling Speed: Delivers rapid and consistent staple deployment.
- Available with a rechargeable battery option to support environmental sustainability.

Ref#	Descripion	Packaging
MNM-EF3-L	E-Force 3 Large	1 PC/BOX
MNM-EF3-M	E-Force 3 Medium	1 PC/BOX
MNM-EF3-S	E-Force 3 Small	1 PC/BOX

MAJnMAR MEDICAL PRODUCTS
203-1130 Austin Ave
Coquitlam, British Columbia
V3K 3P5, CANADA
Tel: +1 778 872 3876
Fax: +1 604 939 9944

Optimal Reloads



Ref #	Size	Number of staples	Staples Height	Description	Packaging
MNM-ECR30GAC	30	48	a=1.0 b= 1.5 c= 2.0	Articulating, 30mm, Grey, with Chip	6 PCS/BOX
MNM-ECR30TAC	30	48	a=2.0 b= 2.5 c= 3.0	Articulating, 30mm, Tan, with Chip	6 PCS/BOX
MNM-ECR30PAC	30	48	a=3.0 b= 3.5 c= 4.0	Articulating, 30mm, Purple, with Chip	6 PCS/BOX
MNM-ECR45GAC	45	66	a=1.0 b= 1.5 c= 2.0	Articulating, 45mm, Grey, with Chip	6 PCS/BOX
MNM-ECR45TAC	45	66	a=2.0 b= 2.5 c= 3.0	Articulating, 45mm, Tan, with Chip	6 PCS/BOX
MNM-ECR45PAC	45	66	a=3.0 b= 3.5 c= 4.0	Articulating, 45mm, Purple, with Chip	6 PCS/BOX
MNM-ECR45BAC	45	66	a=4.0 b= 4.5 c= 5.0	Articulating, 45mm, Black, with Chip	6 PCS/BOX
MNM-ECR60GAC	60	90	a=1.0 b= 1.5 c= 2.0	Articulating, 60mm, Grey, with Chip	6 PCS/BOX
MNM-ECR60TAC	60	90	a=2.0 b= 2.5 c= 3.0	Articulating, 60mm, Tan, with Chip	6 PCS/BOX
MNM-ECR60PAC	60	90	a=3.0 b= 3.5 c= 4.0	Articulating, 60mm, Purple, with Chip	6 PCS/BOX
MNM-ECR60BAC	60	90	a=4.0 b= 4.5 c= 5.0	Articulating, 60mm, Black, with Chip	6 PCS/BOX
MNM-ECR30GATC	30	48	a=1.0 b= 1.5 c= 2.0	Articulating, 30mm, Grey, Tilted Tip, with Chip	6 PCS/BOX
MNM-ECR30TATC	30	48	a=2.0 b= 2.5 c= 3.0	Articulating, 30mm, Tan, Tilted Tip, with Chip	6 PCS/BOX
MNM-ECR30PATC	30	48	a=3.0 b= 3.5 c= 4.0	Articulating, 30mm, Purple, Tilted Tip, with Chip	6 PCS/BOX
MNM-ECR45GATC	45	66	a=1.0 b= 1.5 c= 2.0	Articulating, 45mm, Grey, Tilted Tip, with Chip	6 PCS/BOX
MNM-ECR45TATC	45	66	a=2.0 b= 2.5 c= 3.0	Articulating, 45mm, Tan, Tilted Tip, with Chip	6 PCS/BOX
MNM-ECR45PATC	45	66	a=3.0 b= 3.5 c= 4.0	Articulating, 45mm, Purple, Tilted Tip, with Chip	6 PCS/BOX
MNM-ECR45BATC	45	66	a=4.0 b= 4.5 c= 5.0	Articulating, 45mm, Black, Tilted Tip, with Chip	6 PCS/BOX
MNM-ECR60GATC	60	90	a=1.0 b= 1.5 c= 2.0	Articulating, 60mm, Grey, Tilted Tip, with Chip	6 PCS/BOX
MNM-ECR60TATC	60	90	a=2.0 b= 2.5 c= 3.0	Articulating, 60mm, Tan, Tilted Tip, with Chip	6 PCS/BOX
MNM-ECR60PATC	60	90	a=3.0 b= 3.5 c= 4.0	Articulating, 60mm, Purple, Tilted Tip, with Chip	6 PCS/BOX
MNM-ECR60BATC	60	90	a=4.0 b= 4.5 c= 5.0	Articulating, 60mm, Black, Tilted Tip, with Chip	6 PCS/BOX

MAJnMAR MEDICAL PRODUCTS
203-1130 Austin Ave
Coquitlam, British Columbia
V3K 3P5, CANADA
Tel: +1 778 872 3876
Fax: +1 604 939 9944

PremiaTack Endoscopic Hernia Stapler



FEATURES:

- Our fully disposable Hernia Mesh Fixation Staplers provide Spiral-shaped Titanium Staples and are available in 3 preloaded versions,

5, 20 and 30 staples in 36 cm Length.

- Available in short 23 cm for Pediatric and open surgery applications.

Ref#	Dia	Description	Packaging
MNM-HS15	5mm	Endoscopic Hernia Stapler, 15 Staples, 36cm	1 PC/BOX
MNM-HS20	5mm	Endoscopic Hernia Stapler, 20 Staples, 36cm	1 PC/BOX
MNM-HS30	5mm	Endoscopic Hernia Stapler, 30 Staples, 36cm	1 PC/BOX
MNM-HS15S	5mm	Endoscopic Hernia Stapler, 15 Staples, 23cm	1 PC/BOX
MNM-HS20S	5mm	Endoscopic Hernia Stapler, 20 Staples, 23cm	1 PC/BOX
MNM-HS30S	5mm	Endoscopic Hernia Stapler, 30 Staples, 23cm	1 PC/BOX

PremiaTack A Absorbable Endoscopic Hernia Stapler



FEATURES:

- Absorbable Staple.
- Available with15, 20 or 30 staples in 36 & 23cm Lengths.
- Precise Stapling Mechanism.

Ref#	Dia	Description	Packaging
MNM-HS15A	5mm	Absorbable Hernia Stapler, 15 Staples, 36cm	1 PC/BOX
MNM-HS20A	5mm	Absorbable Hernia Stapler, 20 Staples, 36cm	1 PC/BOX
MNM-HS30A	5mm	Absorbable Hernia Stapler, 30 Staples, 36cm	1 PC/BOX
MNM-HS15AS	5mm	Absorbable Hernia Stapler, 15 Staples, 23cm	1 PC/BOX
MNM-HS20AS	5mm	Absorbable Hernia Stapler, 20 Staples, 23cm	1 PC/BOX
MNM-HS30AS	5mm	Absorbable Hernia Stapler, 30 Staples, 23cm	1 PC/BOX

MAJnMAR MEDICAL PRODUCTS
203-1130 Austin Ave
Coquitlam, British Columbia
V3K 3P5, CANADA
Tel: +1 778 872 3876
Fax: +1 604 939 9944

EndoKnit Suturing Device



FEATURES:

- 360-Degree Cyclic Motion: Features a precision-engineered, slot-driven stitching reload

with smooth 360-degree cyclic motion for consistent performance.

- Angle Adjustment Knob: easily adjusts the horizontal angle of the jaws, making it simple to operate in narrow or restricted spaces.
- Active Release Button: allows for reload release in special situations, reducing surgical risks.
- 60 Degrees Horizontal Articulation: provides adjustable horizontal angles up to 60 degrees for optimal positioning.
- Secure Cap: designed for post-operative cleaning and maintenance.
- Angle Control Lock: ensures precise control over the head and swing positions during use.

Ref#	Description	Length	Packaging
MNM-EK-S	EndoKnit Suturing Device	222mm	1 PC/BOX
MNM-EK-M	EndoKnit Suturing Device	374mm	1 PC/BOX
MNM-EK-L	EndoKnit Suturing Device	449mm	1 PC/BOX
MNM-EKS-2P-200S	Polypropylene EndoKnit Suturing Unit-Self-tied	200mm	1 PC/BOX
MNM-EKS-2P-350S	Polypropylene EndoKnit Suturing Unit-Self-tied	350mm	1 PC/BOX
MNM-EKS-3P-200S	Polypropylene EndoKnit Suturing Unit-Self-tied	200mm	1 PC/BOX
MNM-EKS-3P-350S	Polypropylene EndoKnit Suturing Unit-Self-tied	350mm	1 PC/BOX
MNM-EKS-2P-200	Polypropylene EndoKnit Suturing Unit	200mm	1 PC/BOX
MNM-EKS-2P-350	Polypropylene EndoKnit Suturing Unit	350mm	1 PC/BOX
MNM-EKS-3P-200	Polypropylene EndoKnit Suturing Unit	200mm	1 PC/BOX
MNM-EKS-3P-350	Polypropylene EndoKnit Suturing Unit	350mm	1 PC/BOX
MNM-EKS-2S-200	Silk EndoKnit Suturing Unit	200mm	1 PC/BOX
MNM-EKS-2S-350	Silk EndoKnit Suturing Unit	350mm	1 PC/BOX
MNM-EKS-3S-200	Silk EndoKnit Suturing Unit	200mm	1 PC/BOX
MNM-EKS-3S-350	Silk EndoKnit Suturing Unit	350mm	1 PC/BOX

MAJnMAR MEDICAL PRODUCTS
203-1130 Austin Ave
Coquitlam, British Columbia
V3K 3P5, CANADA
Tel: +1 778 872 3876
Fax: +1 604 939 9944

OptiClip

Disposable Clip Appliers



FEATURES:

- Single-use design ensures sterility and reduces the risk of cross-contamination.
- Ergonomic Handle for comfortable use, reducing surgeon fatigue during prolonged procedures.
- Multiple Clip Sizes: Available in various clip sizes to accommodate different tissue types and surgical needs.

Ref #	Description	Dimension	Packaging
MNM-ECA5	Endoscopic Clip Applier, 5mm, 15 clips	⌀5mmx330mm	1PC/BOX
MNM-ECA5C	Endoscopic Clip Applier, 5mm, 15 clips with Clips counter	⌀5mmx330mm	1PC/BOX
MNM-ECA10	Endoscopic Clip Applier, 10mm, 20 clips, clear shaft	⌀10mmx330mm	1PC/BOX
MNM-ECA10C	Endoscopic Clip Applier, 10mm, 20 clips with Clips counter	⌀10mmx330mm	1PC/BOX

PolyClip Applier

Disposable Polymer Clip Applier + 2 Cartridges:



FEATURES:

- Ergonomic Handle for comfortable use, reducing surgeon fatigue during prolonged procedures.
- Easy Clip Loading: Designed for quick and easy cartridge loading, minimizing downtime during surgery.
- Secure Clip Application: Ensures secure and reliable clip placement for effective ligation or occlusion of vessels and tissues.
- Rotating Shaft: Some models feature a rotating shaft for improved maneuverability and access to hard-to-reach areas.

Ref #	Clip Size	Shaft Diameter (mm)	Working Length	Packaging
MNM-EPCAD-ML	ML	⌀ 5	330	1PC/BOX
MNM-EPCAD-L	L	⌀10	330	1PC/BOX
MNM-EPCAD-XL	XL	⌀10	330	1PC/BOX

MAJnMAR MEDICAL PRODUCTS
203-1130 Austin Ave
Coquitlam, British Columbia
V3K 3P5, CANADA
Tel: +1 778 872 3876
Fax: +1 604 939 9944

PremiaClip

Multi-fire Polymer Clip Appliers



FEATURES:

- Multi-Fire Capability: Designed to apply multiple clips without the need for reloading after each application, improving surgical efficiency.
- Pre-Loaded Clips: allowing for continuous use during surgery.
- A clip counter displaying the remaining number of clips.

Ref #	Clip Size	Shaft Diameter (mm)	Working Length	Clip Quantity	Packaging
MNM-EPCA5-ML	ML	⌀ 5	300	5	1PC/BOX
MNM-EPCA10-ML	ML	⌀ 5	300	10	1PC/BOX
MNM-EPCA15-ML	ML	⌀ 5	300	15	1PC/BOX
MNM-EPCA5-L	L	⌀ 8	290	5	1PC/BOX
MNM-EPCA10-L	L	⌀ 8	290	10	1PC/BOX
MNM-EPCA15-L	L	⌀ 8	290	15	1PC/BOX

FEATURES:

- Multi-Fire Capability allows for the application of multiple clips in a row without needing to remove and reload the instrument.
- Ergonomic Design: designed to reduce hand fatigue during prolonged procedures.
- Lightweight: for easy handle and maneuver.



Ref #	Clip Size	Shaft Diameter (mm)	Working Length	Packaging
MNM-EPCA-15-ML-L	ML	⌀5mm	350mm+/-15mm	1PC/BOX
MNM-EPCA-15-L-L	L	⌀8mm	350mm+/-15mm	1PC/BOX

Reusable Polymer Clip Applier



FEATURES:

- Reusable Design: Made of high-quality, stainless steel that can withstand repeated sterilization.
- Polymer Clip Compatibility: Specifically designed to apply polymer clips.

Ref #	Clip Size	Shaft Diameter (mm)	Working Length
MNM-EPCAR-ML	ML	⌀ 5	330
MNM-EPCAR-L	L	⌀10	330
MNM-EPCAR-XL	XL	⌀10	330

MAJnMAR MEDICAL PRODUCTS
203-1130 Austin Ave
Coquitlam, British Columbia
V3K 3P5, CANADA
Tel: +1 778 872 3876
Fax: +1 604 939 9944

PolyClip

Polymer Clip



FEATURES:

- Biocompatibility: Made of medical-grade polymers that are safe for use in the human body.
- Radiolucency: Unlike metal clips, polymer clips are radiolucent, meaning they do not interfere with imaging techniques such as X-rays, CT scans,

or MRIs.

- Ideal for patients who may require postoperative imaging.
- Lightweight: Polymer clips are lightweight, making them easier to handle and less likely to cause tissue damage.

Ref #	Clip Size	Length	Color	Packaging
MNM-PC-ML	ML	9.5	Green	15Cartridges/6 per Box
MNM-PC-L	L	13	Purple	15Cartridges/6 per Box
MNM-PC-XL	XL	17	Yellow	15Cartridges/6 per Box

Reusable Clip Remover

FEATURES:



- Reusable Design: Designed for multiple uses, making it a cost-effective option for surgical procedures.
- Made from high-quality, medical-grade materials to ensure durability and reliability.
- Ergonomic Handle: Features an ergonomic handle for comfortable use, reducing surgeon fatigue during prolonged

procedures.

Ref #	Clip Size	Shaft Diameter (mm)	Working Length
MNM-PCR-5	ML	⌀ 5	330
MNM-PCR-10	L/XL	⌀ 10	330

MAJnMAR MEDICAL PRODUCTS
203-1130 Austin Ave
Coquitlam, British Columbia
V3K 3P5, CANADA
Tel: +1 778 872 3876
Fax: +1 604 939 9944

Reusable Titanium Clip Appliers

FEATURES:

- Secure Clip Application: Ensures secure and reliable clip placement for effective ligation or occlusion of vessels and tissues.
- Laparoscopic and Open Surgery Compatibility: Suitable for both laparoscopic (minimally invasive) and open surgical procedures.
- High-Quality Materials: Made of high-quality, medical-grade materials to ensure durability and reliability.

	Ref #	Description	Dimension	Packaging
	MNM-221-05116	Small clip applier, single-action jaw	⌀10mmx330mm	1PC/BOX
	MNM-221-05126	Medium clip applier, single-action jaw	⌀5mmx330mm	1PC/BOX
	MNM-221-05136	Medium-Large clip applier, single-action jaw	⌀10mmx330mm	1PC/BOX
	MNM-221-05146	Large clip applier applier, single- action jaw	⌀10mmx330mm	1PC/BOX
	MNM-221-00516	Small clip applier for open surgery	18mm/28mm	1PC/BOX
	MNM-221-00526	Medium clip applier for open surgery	18mm/28mm	1PC/BOX
	MNM-221-00536	Medium-large clip applier for open surgery	23mm/28mm	1PC/BOX
	MNM-221-00546	Large clip applier for open surgery	28mm	1PC/BOX

Titanium Clips

FEATURES:

- A diamond-shaped inner clip profile provides an enlarged contact surface.
- The special closing characteristic of the clip allows a repositioning of the clip and reduces the risk of tissue slippage.

	Ref#	Description	Dimension	Packaging
	MNM-TC-S	Small Titanium Clip	S	36 Cartridges/BOX
	MNM-TC-M	Medium Titanium Clip	M	36 Cartridges/BOX
	MNM-TC-ML	Medium Large Titanium Clip	M/L	20 Cartridges/BOX
	MNM-TC-L	Large Titanium Clip	L	20 Cartridges/BOX

MAJnMAR MEDICAL PRODUCTS
203-1130 Austin Ave
Coquitlam, British Columbia
V3K 3P5, CANADA
Tel: +1 778 872 3876
Fax: +1 604 939 9944

Veress Needle



FEATURES:

- Sharp bevelled tip ensures quick, easy and safe controlled insertion.
- Red safety indicator.
- Spring-loaded style provides audible "click" when piercing tissue layers.

Ref#	Dimension	Packaging
MNM-121-20127	2mmx120mm	5PCS/BOX
MNM-121-20157	2mmx150mm	5PCS/BOX

Suction irrigation



FEATURES:

- Ergonomic design for easy of use.
- One way valve prevents the saline from refluxing.
- Differentiated color buttons and tube

for easy identification.

- Smooth movement of the handpiece prevents fatigue during the operation.

Ref#	Description	Dimension	Packaging
MNM-121-33004	Suction & Irrigation Tube	⌀5x330	5PCS/BOX

Endobag



FEATURES:

- Specimen Retrieval Bag.
- Unique ergonomic design offers ease of use.
- 10mm introducer provides ease of insertion.
- Self-opening design provides easy and safe specimen retrieval.
- High strength thermoplastic polyurethane pouch material prevents leakage.
- Enhanced titanium ring ensures large opening and good fitting in specimen loading.

Ref#	Description	Dimension	Packaging
MNM-121-00318	Specimen Retrieval Bag without string 100ml	100ml	20PCS/BOX
MNM-121-00328	Specimen Retrieval Bag without string 200m	200ml	20PCS/BOX
MNM-121-00348	Specimen Retrieval Bag without string 400ml	400ml	20PCS/BOX
MNM-121-00378	Specimen Retrieval Bag without string 700ml	700ml	20PCS/BOX
MNM-121-01318	Specimen Retrieval Bag with string 100ml	100ml	20PCS/BOX
MNM-121-01328	Specimen Retrieval Bag with string 200ml	200ml	20PCS/BOX
MNM-121-01348	Specimen Retrieval Bag with string 400ml	400ml	20PCS/BOX
MNM-121-01378	Specimen Retrieval Bag with string 700ml	700ml	20PCS/BOX

MAJnMAR MEDICAL PRODUCTS
203-1130 Austin Ave
Coquitlam, British Columbia
V3K 3P5, CANADA
Tel: +1 778 872 3876
Fax: +1 604 939 9944

Needle Holder

FEATURES:

- Ultra-light design greatly reduce the fatigue for long life of use.
- Excellent titanium material and precise spring contribute long life cycle.
- Variable tips and handles selection .

V Handle



Ref#	Description	Dimension	Packaging
MNM-221-56616	V Handle Straight Needle Holder	⌀5x330mm	1PC/BOX
MNM-221-56626	V Handle Left Curved Needle Holder	⌀5x330mm	1PC/BOX
MNM-221-56636	V Handle Right Curved Needle Holder	⌀5x330mm	1PC/BOX
MNM-221-56646	V Handle Reposition Needle Holder	⌀5x330mm	1PC/BOX
MNM-221-56656	V Handle Multifunctional Reposition Needle Holder	⌀5x330mm	1PC/BOX
MNM-221-56716	Gun-Shaped Handle Straight Needle Holder	⌀5x330mm	1PC/BOX
MNM-221-56726	Gun-Shaped Handle Left Curved Needle Holder	⌀5x330mm	1PC/BOX
MNM-221-56736	Gun-Shaped Handle Right Curved Needle Holder	⌀5x330mm	1PC/BOX
MNM-221-56746	Gun-Shaped Handle Needle Holder with Reposition Function	⌀5x330mm	1PC/BOX
MNM-221-56756	Gun-Shaped Handle Multifunctional Reposition Needle Holder	⌀5x330mm	1PC/BOX

Gun Handle



Ref#	Description	Dimension	Packaging
MNM-221-56716	Gun-Shaped Handle Straight Needle Holder	⌀5x330mm	1PC/BOX
MNM-221-56726	Gun-Shaped Handle Left Curved Needle Holder	⌀5x330mm	1PC/BOX
MNM-221-56736	Gun-Shaped Handle Right Curved Needle Holder	⌀5x330mm	1PC/BOX
MNM-221-56746	Gun-Shaped Handle Needle Holder with Reposition Function	⌀5x330mm	1PC/BOX
MNM-221-56756	Gun-Shaped Handle Multifunctional Reposition Needle Holder	⌀5x330mm	1PC/BOX

Endoscopic Instruments



MAJnMAR MEDICAL PRODUCTS
203-1130 Austin Ave
Coquitlam, British Columbia
V3K 3P5, CANADA
Tel: +1 778 872 3876
Fax: +1 604 939 9944

Disposable Instruments

FEATURES:


- A broad range of Laparoscopic Instruments for wide surgical applications.

MNM-302.802 
MNM-302.812 
MNM-302.822 
MNM-302.823 
MNM-302.839 

	Code	Description	Dimension	Pack Qty
	MNM-302.811	Straight Metzenbaum Scissors	⌀5mm x 330mm	12Pcs/Box
	MNM-302.812	Curved Metzenbaum Scissors	⌀5x330mm , blade length 16.5mm	12pcs/box
	MNM-302.813	Curved Metzenbaum Micro Scissors	⌀5x330mm , blade length 13.5mm	12pcs/box
	MNM-302.815	Fully Insulated Curved Metzenbaum Scissors	⌀5x330mm , blade length 15.5mm	12pcs/box
	MNM-302.816	Fully Insulated Curved Micro Metzenbaum Scissors	⌀5x330mm , blade length 10mm	12pcs/box
	MNM-302.802	Curved Maryland Dissecting Forceps	⌀5x330mm	12pcs/box
	MNM-302.805	Modified Maryland Dissector (90°)	⌀5x330mm	12pcs/box
	MNM-302.821	Allis Grasping Forceps	⌀5x330mm	12pcs/box
	MNM-302.822	Fenestrated Johan Grasping Forceps	⌀5x330mm	12pcs/box
	MNM-302.823	Clinch Grasping Forceps	⌀5x330mm	12pcs/box
	MNM-302.838	Single Action Fenestrated Johan Grasping Forceps	⌀5x330mm jaw length 22mm	12pcs/box
	MNM-302.839	Single Action Fenestrated Johan Grasping Forceps	⌀5x330mm jaw length 40mm	12pcs/box
	MNM-302.840	Wave Teeth Grasping Forceps	⌀5x330mm	12pcs/box
	MNM-302.841	Atraumatic Grasping Forceps	⌀5x330mm	12pcs/box
	MNM-302.842	Babcock Grasping Forceps	⌀5x330mm	12pcs/box
	MNM-302.945	Claw Grasping Forceps	⌀5x330mm	12pcs/box
	MNM-302.946	Duck Jaw Grasping Forceps	⌀5x330mm	12pcs/box
	MNM-302.901	Bipolar Single Action Maryland Dissecting Forceps	⌀5x330mm	10pcs/box
	MNM-121-30213	Cholangiography Forceps	⌀5x330mm	12pcs/box

MAJnMAR MEDICAL PRODUCTS
203-1130 Austin Ave
Coquitlam, British Columbia
V3K 3P5, CANADA
Tel: +1 778 872 3876
Fax: +1 604 939 9944

Coagulation Electrode



FEATURES:

- Convenient, single-use instrument guarantees sterility and fasten monopolar electrosurgical activation.
- Fully insulated 330mm shaft provides safe current containment .
- Lightweight, ergonomic handle.

	Ref #	Description	Packaging
	MNM-122-30220	Ball, ⌀5mm x 330mm	30PCS/BOX
	MNM-122-30110	J-Hook, 5mm x 330mm	30PCS/BOX
	MNM-122-30440	Spatula, ⌀5mmx330mm	30PCS/BOX
	MNM-122-30330	Needle, ⌀5mmx330mm	30PCS/BOX
	MNM-122-30660	L-Hook, ⌀5mmx330mm	30PCS/BOX

Blunt Tip Dissector, Knot Pusher, Camera Sleeve

FEATURES:

- Blunt Tip Dissector: For safe tissue dissection and manipulation.
- Knot Pusher: For securing and tightening sutures during minimally invasive surgeries.
- Camera Sleeve: For protecting the laparoscopic camera and maintaining sterility.

	Ref #	Description	Demnision	Packaging
	MNM-121-S4559	BluntTip Dissector	⌀5mmx400m	30PCS/BOX
	MNM-21-30119	Knot Pusher	⌀5mmx330mm	30PCS/BOX
	MNM-121-00449	Camera Sleeve	Length 3000m	25PCS/BOX

MAJnMAR_MEDICAL_PRODUCTS
203-1130 Austin Ave
Coquitlam, British Columbia
V3K 3P5, CANADA
Tel: +1 778 872 3876
Fax: +1 604 939 9944

Linear Cutting Staplers & Reloads




FEATURES:

- Simultaneous Cutting and Stapling: Combines the functions of cutting and stapling in a single device, improving surgical efficiency.
- Offers variable staple heights to accommodate varying tissue thicknesses.

Ref#	Description	Packaging
MNM-LC60W	60mm cutting length, 2.5mm staple height, White	1PC/BOX
MNM-LC60B	60mm cutting length, 3.8mm staple height, Blue	1PCS/BOX
MNM-LC60G	60mm cutting length, 4.8mm staple height, Green	1PC/BOX
MNM-LC80B	80mm cutting length, 3.8mm staple height, Blue	1PC/BOX
MNM-LC80D	80mm cutting length, 4.2mm staple height, Gold	1PC/BOX
MNM-LC80G	80mm cutting length, 4.8mm staple height, Green	1PC/BOX
MNM-LC100B	100mm cutting length, 3.8mm staple height, Blue	1PC/BOX
MNM-LC100G	100mm cutting length, 4.8mm staple height, Green	1PC/BOX

FEATURES:

- Enhanced Hemostasis: Provides excellent hemostatic control, minimizing bleeding during and after stapling.
- Ease of Loading: Designed for quick and easy loading into the Linear cutting Stapler, minimizing downtime during surgery.



Ref#	Description	Packaging
MNM-LCR60W	60mm cutting length, 2.5mm staple height, White	6PCS/BOX
MNM-LCR60B	60mm cutting length, 3.8mm staple height, Blue	6PCS/BOX
MNM-LCR60G	60mm cutting length, 4.8mm staple height, Green	6PCS/BOX
MNM-LCR80B	80mm cutting length, 3.8mm staple height, Blue	6PCS/BOX
MNM-LCR80D	80mm cutting length, 4.2mm staple height, Gold	6PCS/BOX
MNM-LCR80G	80mm cutting length, 4.8mm staple height, Green	6PCS/BOX
MNM-LCR100B	100mm cutting length, 3.8mm staple height, Blue	6PCS/BOX
MNM-LCR100G	100mm cutting length, 4.8mm staple height, Green	6PCS/BOX

MAJnMAR MEDICAL PRODUCTS
203-1130 Austin Ave
Coquitlam, British Columbia
V3K 3P5, CANADA
Tel: +1 778 872 3876
Fax: +1 604 939 9944

Linear Staplers & Reloads

	Ref #	Description	Packaging
	MNM-LS30B	30mm stapling length, 3.8mm staple height, Blue	1PC/BOX
	MNM-LS30G	30mm stapling length, 4.8mm staple height, Green	1PC/BOX
	MNM-LS45B	45mm stapling length, 3.8mm staple height, Blue	1PC/BOX
	MNM-LS45G	45mm stapling length, 4.8mm staple height, Green	1PC/BOX
	MNM-LS60B	60mm stapling length, 3.8mm staple height, Blue	1PC/BOX
	MNM-LS60G	60mm stapling length, 4.8mm staple height, Green	1PC/BOX
	MNM-LSR30B	30mm stapling length, 3.5mm staple height, Blue	6PCS/BOX
	MNM-LSR30G	30mm stapling length, 4.8mm staple height, Green	6PCS/BOX
	MNM-LSR45B	45mm stapling length, 3.5mm staple height, Blue	6PCS/BOX
	MNM-LSR45G	45mm stapling length, 4.8mm staple height, Green	6PCS/BOX
	MNM-LSR60B	60mm stapling length, 3.5mm staple height, Blue	6PCS/BOX
	MNM-LSR60G	60mm stapling length, 4.8mm staple height, Green	6PCS/BOX


MAJnMAR MEDICAL PRODUCTS
203-1130 Austin Ave
Coquitlam, British Columbia
V3K 3P5, CANADA
Tel: +1 778 872 3876
Fax: +1 604 939 9944

Circular Stapler



INDICATIONS:

- The circular stapler is used to perform end to side, side to end and end to end anastomosis canal organs such as esophagus, stomach, and intestine.



Ref #	MNM-CS24-3	MNM-CS26-3	MNM-CS29-3	MNM-CS32-3
Number of rows	3	3	3	3
Stapler Material	Titanium - TA1G	Titanium - TA1G	Titanium - TA1G	Titanium - TA1G
Total length of the stapler (mm)	448±10	448±10	448±10	448±10
External diameter of the staples cartridge sleeve (mm)	24.1±2	26.1±2	29.1±2	32.1±2
External diameter of the cutting knife (mm)	15±1	17±1	20.3±1	22±1
Number of Staples	30	36	36	45
Staples height (mm)	4.8mm	4.8mm	4.8mm	5.0mm
Staples wire diameter (mm)	0.3mm	0.3mm	0.3mm	0.3mm
Outer diameter of the guide bushing (mm)	6.12mm	6.12mm	6.12mm	6.12mm
Total length of the guide bushing (mm)	31.84mm	31.84mm	31.84mm	31.84mm
Outer diameter of the extension (mm)	6.3mm	6.3mm	6.3mm	6.3mm
Total length of the extension (mm)	49.76mm	49.76mm	49.76mm	49.76mm
Outer diameter of the extension (mm)	7.9mm	7.9mm	7.9mm	7.9mm
Total length of the extension (mm)	118.82mm	118.82mm	118.82mm	118.82mm
Color	White	Gold	Purple	Black

Ref #	MNM-CS21	MNM-CS25	MNM-CS29	MNM-CS32	MNM-CS34
Number of rows	2	2	2	2	2
Staples Material	Titanium - TA1G	Titanium - TA1G	Titanium - TA1G	Titanium - TA1G	Titanium - TA1G
The total length of the stapler (mm).	445±10	445±10	445±10	445±10	445±10
External diameter of the staples cartridge sleeve (mm)	22.1±2	26.3±2	29±2	32±2	33.9±2
Diameter of the nail seat	21.6±1	25.6±1	28.6±1	31.6±1	33.6±1
External diameter of the cutting knife (mm)	13±1	17±1	20±1	22±1	24±1
Number of Staples	18	20	24	30	32
Stapler height (mm)	4.5mm	4.8mm	4.8mm	5.0mm	5.1mm
Staple diameter (mm)	0.28±0.2mm	0.28±0.2mm	0.28±0.2mm	0.28±0.2mm	0.28±0.2mm
Outer diameter of the guide bushing(mm)			6.1±0.5mm		
Total length of the guide bushing (mm)			31.7±3mm		
Outer diameter of the extension part (mm)			6.2±0.5mm		
Total length of the extension part (mm)			49.8±5mm		

MAJnMAR MEDICAL PRODUCTS
203-1130 Austin Ave
Coquitlam, British Columbia
V3K 3P5, CANADA
Tel: +1 778 872 3876
Fax: +1 604 939 9944

Ref #	MNM-CS21	MNM-CS25	MNM-CS29	MNM-CS32	MNM-CS34
Outer diameter of the extension piece (mm)			8±2mm		
Total length of the extension piece (mm)			115.2±5mm		
Color	Dark blue	Light blue	Green	White	Grey

Hemorrhoids Stapler




INDICATIONS:

- Hemorrhoids Stapler is used for excision of mucosal and submucosal tissues along the dentate line, to restore the normal anatomical structure of lower rectum.

Ref #	MNM-HS32	MNM-HS34
Number of rows	2	2
Staples Material	Titanium-TA1G	Titanium-TA1G
The total length of the stapler (in mm).	445±10	445±10
The external diameter of the staples cartridge sleeve (mm)	32±2	34±2
External diameter of the cutting knife (mm)	23±1	24±1
Number of Staples	32	34
Staple height	4.0±0.2	4.0±0.2
Staples wire diameter (mm)	0.3mm	0.3mm
Color	White	Green
Closed Staple Shape	B	B
Closed Staples Height	0.75-1.5mm	0.75-1.5mm
Ref #	MNM-HS32-3	MNM-HS34-3
Number of rows	3	3
Staples Material	Titanium-TA2G	Titanium-TA2G
Total length of the stapler (mm)	386±10	386±10
The external diameter of the staples cartridge sleeve (mm)	32.9±1.0	34.5±1.0
External diameter of the cutting knife (mm)	22.6±0.5	23.4±0.5
Number of Staples	48	48
Staples height (mm)	3.8±0.3	3.8±0.3
Staples wire diameter (mm)	0.3mm	0.3mm
Color	White	Green

MAJnMAR MEDICAL PRODUCTS
203-1130 Austin Ave
Coquitlam, British Columbia
V3K 3P5, CANADA
Tel: +1 778 872 3876
Fax: +1 604 939 9944

Silicone Drainage System



Made of 100% medical grade silicone, we produce a complete line of Silicone Closed Wound Drains and Reservoirs for multiple surgical applications.

Silicone provides superior biocompatibility, flexibility, and great comfort to patients. MAJnMAR drains deliver extremely low tissue responses to wounds, minimizing tissue-adhesion and tissue-trauma during drain removal.

MAJnMAR MEDICAL PRODUCTS
203-1130 Austin Ave
Coquitlam, British Columbia
V3K 3P5, CANADA
Tel: +1 778 872 3876
Fax: +1 604 939 9944

Round channeled Drain



FEATURES:

- Multiple narrow ducts structured for maximized drainage.
- Reduce withdrawal difficulty against tissue ingrowth.
- Available in four sizes (10, 15, 19 & 24 Fr).

	Type	Ref #	Description	Packaging
	Set	MNM-RDF10T	Drain 10Fr with Trocar	10PCS/BOX
		MNM-RDF15T	Drain 15Fr with Trocar	10PCS/BOX
		MNM-RDF19T	Drain 19Fr with Trocar	10PCS/BOX
		MNM-RDF24T	Drain 24Fr with 16Fr Tube	10PCS/BOX
		MNM-RDF24TT	Drain 24Fr with 16Fr Tube and Trocar	10PCS/BOX
	Kit	MNM-RDF10T-RXX	Drain 10Fr, Trocar & Reservoir 100/150/200/400cc	10PCS/BOX
		MNM-RDF15T-RXX	Drain 15Fr, Trocar & Reservoir 100/150/200/400cc	10PCS/BOX
		MNM-RDF19T-RXX	Drain 19Fr, Trocar & Reservoir 100/150/200/400cc	10PCS/BOX
		MNM-RDF24T-RXX	Drain 24Fr, Trocar & Reservoir 100/150/200/400cc	10PCS/BOX
	System	RDF10T-RXX-BYY	Drain 10Fr, Trocar, Reservoir 100/150/200/400cc & Bag 1000/2000ml	10PCS/BOX
		RDF15T-RXX-BYY	Drain 15Fr, Trocar, Reservoir 100/150/200/400cc & Bag 1000/2000ml	10PCS/BOX
		RDF19T-RXX-BYY	Drain 19Fr, Trocar, Reservoir 100/150/200/400cc & Bag 1000/2000ml	10PCS/BOX
		RDF24T-RXX-BYY	Drain 24Fr, Trocar, Reservoir 100/150/200/400cc & Bag 1000/2000ml	10PCS/BOX

MAJnMAR MEDICAL PRODUCTS
203-1130 Austin Ave
Coquitlam, British Columbia
V3K 3P5, CANADA
Tel: +1 778 872 3876
Fax: +1 604 939 9944

Round Perforated Drain



FEATURES:

- Multiple perforations for maximized drainage.
- Available in 7 sizes (8, 10, 12, 14, 16, 18 & 20 Fr).

	Type	Ref #	Description	Packaging
	Set	MNM-RDP08T	Drain 8Fr with Trocar	10PCS/BOX
		MNM-RDP10T	Drain 10Fr with Trocar	10PCS/BOX
		MNM-RDP12T	Drain 12Fr with Trocar	10PCS/BOX
		MNM-RDP14T	Drain 14Fr with Trocar	10PCS/BOX
		MNM-RDP16T	Drain 16Fr with Trocar	10PCS/BOX
		MNM-RDP18T	Drain 18Fr with Trocar	10PCS/BOX
		MNM-RDP20T	Drain 20Fr with Trocar	10PCS/BOX
	Kit	MNM-RDP08T-RXX	Drain 8Fr, Trocar & Reservoir 100/150/200/400cc	10PCS/BOX
		MNM-RDP10T-RXX	Drain 10Fr, Trocar & Reservoir 100/150/200/400cc	10PCS/BOX
		MNM-RDP12T-RXX	Drain 12Fr, Trocar & Reservoir 100/150/200/400cc	10PCS/BOX
		MNM-RDP14T-RXX	Drain 14Fr, Trocar & Reservoir 100/150/200/400cc	10PCS/BOX
		MNM-RDP16T-RXX	Drain 16Fr, Trocar & Reservoir 100/150/200/400cc	10PCS/BOX
		MNM-RDP18T-RXX	Drain 18Fr, Trocar & Reservoir 100/150/200/400cc	10PCS/BOX
		MNM-RDP20T-RXX	Drain 20Fr, Trocar & Reservoir 100/150/200/400cc	10PCS/BOX
	System	RDP08T-RXX-BYY	Drain 8Fr, Trocar, Reservoir 100/150/200/400cc & Bag 1000/2000 ml	10PCS/BOX
		RDP10T-RXX-BYY	Drain 10Fr, Trocar, Reservoir 100/150/200/400cc & Bag 1000/2000 ml	10PCS/BOX
		RDP12T-RXX-BYY	Drain 12Fr, Trocar, Reservoir 100/150/200/400cc & Bag 1000/2000 ml	10PCS/BOX
		RDP14T-RXX-BYY	Drain 14Fr, Trocar, Reservoir 100/150/200/400cc & Bag 1000/2000 ml	10PCS/BOX
		RDP16T-RXX-BYY	Drain 16Fr, Trocar, Reservoir 100/150/200/400cc & Bag 1000/2000 ml	10PCS/BOX
		RDP18T-RXX-BYY	Drain 18Fr, Trocar, Reservoir 100/150/200/400cc & Bag 1000/2000 ml	10PCS/BOX
		RDP20T-RXX-BYY	Drain 20Fr, Trocar, Reservoir 100/150/200/400cc & Bag 1000/2000 ml	10PCS/BOX

Flat Channeled Drain



FEATURES:

- Multiple narrow ducts structure for maximized drainage.
- Available in 2 sizes (8, 11mm).

	Type	Ref #	Description	Packaging
	Set	MNM-FDF08T	Drain 8mmwith Trocar	10PCS/BOX
		MNM-FDF11T	Drain 11mm with Trocar	10PCS/BOX
	Kit	MNM-FDF08T-RXX	Drain 8mm, Trocar & Reservoir 100/150/200/400cc	10PCS/BOX
		MNM-FDF11T-RXX	Drain 11mm, Trocar & Reservoir 100/150/200/400cc	10PCS/BOX
	System	FDF08T-RXX-BYY	Drain 8mm, Trocar & Reservoir 100/150/200/400cc & Bag 800/1000/2000ml	10PCS/BOX
		FDF11T-RXX-BYY	Drain 11mm, Trocar & Reservoir 100/150/200/400cc & Bag 800/1000/2000ml	10PCS/BOX

MAJnMAR MEDICAL PRODUCTS
203-1130 Austin Ave
Coquitlam, British Columbia
V3K 3P5, CANADA
Tel: +1 778 872 3876
Fax: +1 604 939 9944

Flat Perforated Drain



FEATURES:

- Multiple perforations for maximized drainage.
- Available in 3 sizes (7, 10 & 13 mm).

	Type	Ref #	Description	Packaging
	Set	MNM-FDP07T	Drain 7mm with Trocar	10PCS/BOX
		MNM-FDP10T	Drain 10mm with Trocar	10PCS/BOX
		MNM-FDP13T	Drain 13mm with Trocar,	10PCS/BOX
	Kit	MNM-FDP07T-RXX	Drain 7mm , Trocar & Reservoir 100/150/200/400cc	10PCS/BOX
		MNM-FDP10T-RXX	Drain 10mm, Trocar & Reservoir 100/150/200/400cc	10PCS/BOX
		MNM-FDP13T-RXX	Drain 13mm , Trocar & Reservoir 100/150/200/400cc	10PCS/BOX
	System	FDP07T-RXX-BYY	Drain 7mm, Trocar & Reservoir 100/150/200/400cc & Bag 800/1000/2000ml	10PCS/BOX
		FDP10T-RXX-BYY	Drain 10mm, Trocar & Reservoir 100/150/200/400cc & Bag 800/1000/2000ml	10PCS/BOX
		FDP13T-RXX-BYY	Drain 13mm, Trocar & Reservoir 100/150/200/400cc & Bag 800/1000/2000ml	10PCS/BOX

Silicone Bulb reservoir



FEATURES:

- Light and patient-friendly system used with closed wound drainage offers a constant suction profile.
- Resulting in improved tissue approximation and reduced risk of hematoma.

	Ref #	Description	Packaging
	MNM-R100	Bulb Reservoir 100CC	10PCS/BOX
	MNM-R150	Bulb Reservoir 150CC	8PCS/BOX
	MNM-R200	Bulb Reservoir 200CC	8PCS/BOX
	MNM-R400	Bulb Reservoir 400CC	10PCS/BOX
	MNM-R100B	Bulb Reservoir with Bottom Outlet 100CC	10PCS/BOX
	MNM-R150B	Bulb Reservoir with Bottom Outlet 150CC	8PCS/BOX
	MNM-R200B	Bulb Reservoir with Bottom Outlet 200CC	8PCS/BOX
	MNM-R400B	Bulb Reservoir with Bottom Outlet 400CC	10PCS/BOX

MAJnMAR MEDICAL PRODUCTS
203-1130 Austin Ave
Coquitlam, British Columbia
V3K 3P5, CANADA
Tel: +1 778 872 3876
Fax: +1 604 939 9944

I -Vac Reservoir



FEATURES:

- A sterile, disposable, and portable system used with closed wound drainage that offers a constant suction profile resulting in improved tissue approximation and reduced risk of hematoma.
- The system features a pre-compressed state with a standard Y-connector for dual drainage.

Ref#	Description	Packaging
MNM-S23102	i-Vac Reservoir 300mL	10PCS/BOX
MNM-S23103	i-Vac Reservoir 450mL	10PCS/BOX

Silicone T-Tube Drain



FEATURES:

- Material: Made of high-quality medical-grade silicone, which is biocompatible, flexible, and resistant to kinking.
- T-Shaped Design: Features a T-shaped structure with a long vertical limb and two shorter horizontal limbs, allowing it to anchor securely within a duct or tube.
- Soft and Flexible: silicone material, reducing the risk of tissue trauma and patient discomfort.

Ref#	Description	Packaging
MNM-TT10	Silicone T-Tube Drain 10Fr	10PCS/BOX
MNM-TT12	Silicone T-Tube Drain 12Fr	10PCS/BOX
MNM-TT14	Silicone T-Tube Drain 14Fr	10PCS/BOX

MAJnMAR MEDICAL PRODUCTS
203-1130 Austin Ave
Coquitlam, British Columbia
V3K 3P5, CANADA
Tel: +1 778 872 3876
Fax: +1 604 939 9944

Calibration Tube

CT1

Consists of a catheter, an internal conical connector, a silicone head a universal connector and a flow stop plate.



Ref #	Size (Fr)	Length (mm)	ID (mm)	OD (mm)
MNM-CT132	32	1200	6.7	10.7
MNM-CT134	34	1200	7.1	11.3
MNM-CT136	36	1200	8.0	12
MNM-CT138	38	1200	8.4	12.4
MNM-CT140	40	1200	8.2	13.2

CT2

Consists of a catheter, an internal conical connector, a silicone head a universal connector, a flow stop plate, a one-way valve and a balloon.



Ref #	Size (Fr)	Length (mm)	ID (mm)	OD (mm)	Balloon Dia.	Balloon Vol.
MNM-CT232	32	745	6.7	10.7	41	30
MNM-CT234	34	745	7.1	11.3	41	30
MNM-CT236	36	745	8.0	12	41	30
MNM-CT238	38	745	8.4	12.4	41	30
MNM-CT240	40	745	8.2	13.2	41	30